

## **IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY**

**DANGER-** TO REDUCE THE RISK OF BURNS, ELECTRIC SHOCK AND FIRE, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS.

1. READ ALL INSTRUCTIONS BEFORE USING.
2. DO NOT USE WHILE SLEEPING. DO NOT USE ON AN INFANT.
3. THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, A PARALYZED PERSON, A SLEEPING OR

UNCONSCIOUS PERSON, A PERSON WITH DIABETES, OR A PERSON WITH POOR BLOOD CIRCULATION. DO NOT USE ON AREAS OF INSENSITIVE SKIN.

4. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING. CHECK SKIN UNDER PAD FREQUENTLY TO AVOID BURNING AND BLISTERING.
5. PLACE PAD ON TOP OF AND NOT UNDER THE PART OF THE BODY NEEDING HEAT - BURNS CAN RESULT REGARDLESS OF CONTROL SETTING.
6. DO NOT USE IN OXYGEN EN-

RICHED ENVIRONMENT OR NEAR EQUIPMENT WHICH STORES OR EMITS OXYGEN.

7. **DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.**
8. **DO NOT SIT ON OR AGAINST OR CRUSH PAD - AVOID SHARP FOLDS.**
9. **THIS POWER UNIT IS INTENDED TO BE CORRECTLY ORIENTATED IN A VERTICAL OR FLOOR MOUNT POSITION.**
10. **NEVER PULL THIS PAD BY THE SUPPLY CORD AND DO NOT USE THE CORD AS A HANDLE.**

11. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS.

12. DISCARD THE PAD IF COVER SHOWS ANY SIGNS OF DETERIORATION (SUCH AS BLISTERING OR CRACKING).

13. USE THE SUPPLY TRANSFORMER ONLY ON A 120 VOLT AC CIRCUIT. UNPLUG WHEN NOT IN USE.

14. DO NOT TAMPER WITH THIS PAD IN ANY WAY. THERE ARE NO USER SERVICEABLE PARTS. IF FOR ANY REASON THIS PAD DOES NOT FUNCTION SATISFAC-

TORILY, RETURN TO AUTHORIZED SERVICE CENTER AT ADDRESS GIVEN.

15. DO NOT USE THIS PAD WITH LINIMENT, SALVE, OR OINTMENT PREPARATIONS THAT CONTAIN HEAT PRODUCING INGREDIENTS. SKIN BURNS COULD RESULT.

16. USE ONLY SUPPLIED TRANSFORMER WITH THIS PAD.

17. SAVE THESE INSTRUCTIONS.

## CAUTIONS

- This massager should not be used in the following cases unless directed by a physician: Tuberculosis, benign and malignant tumors, hemorrhages, inflammations of the skin, phlebitis and thrombosis, open and fresh wounds, bruises, broken skin or varicose veins or undetermined calf pain or joint pain.
- Consult your physician, should you have any therapeutic questions regarding its use.
- In case of pregnancy, diabetes, or illness, consult your physician before using massager.
- Do not use on any unexplained pain before consulting your physician.
- Massage should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your physician.
- Do not fall asleep while using this massager.
- Keep away from children. This unit is not a toy. Children should not use it or play with it.

## MAINTENANCE INSTRUCTIONS

- 1) Your massager is designed for household use, and is virtually maintenance-free. No lubrication is needed. To clean the head and body of the massager, use a damp cloth only. Do not use any strong chemical cleaners such as Benzine.
- 2) If the cord or massager is damaged or does not operate properly, return it to an authorized service center. Any required servicing should be performed by a Conair authorized service centre.
- 3) Never allow the power supply cord to be excessively pulled or twisted. Inspect the cord frequently for damage. Stop use immediately if damage is visible, or unit stops, or operates intermittently. Never wrap the cord around the appliance.

## RECOMMENDED MASSAGING TIMES

- Generally, massage should last 10 - 20 minutes on any muscle of the body. Treatment may be applied once or twice daily.
- As people's metabolisms varies, so does their response to massage and you will soon adjust the time to suit your particular needs and personal comfort. If you have any doubts regarding treatment, we recommend you consult your physician.

## TO OPERATE

To operate the massaging heating pad, first set both the massaging and heat switches to off and plug the transformer into a standard 120v electrical outlet. Plug the power jack into the hand controller. Now turn on either the massaging feature by switching it to HI or LO or the heat feature to "ON". You can use these 2 features together or alone.

For maximum benefit, use should not exceed 20 minutes at one time.

## TIMER FEATURE

The massage motors will automatically shut-off after approximately 20 minutes. To restart, first set the massage switch back to OFF and then set the switch to LO or HI again. Please note the heating pad operates independently of the timer.

## GET TO KNOW YOUR MASSAGING HEATING PAD

King Size Pad  
11.5" x 24"

Self-adhesive Belt  
for easy application

CUL Listed  
Transformer  
Included

On/Off Heat and  
Off/Hi/Low Massage  
Option

Massaging Heating Pad  
for custom comfort to ease  
muscle aches and pains



## ONE YEAR LIMITED WARRANTY

Conair will (at our option) repair or replace your appliance free of charge for 12 months from the date of purchase if the appliance is defective in workmanship or materials. To obtain service under this warranty, return the defective product to the service center nearest you, together with your sales slip and \$7.50 for postage and handling. In the absence of a purchase receipt, the warranty period will be 12 months from the date of manufacture.

**ANY IMPLIED WARRANTIES, OBLIGATIONS OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 12 MONTH DURATION OF THIS WARRANTY.**

**IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER.**

This warranty gives you specific legal rights, and you may also have other rights which vary from province to province.

## CONAIR CONSUMER PRODUCTS ULC

### SERVICE CENTRE:

Conair Consumer Products ULC  
100 Conair Parkway,  
Woodbridge, ON L4H 0L2

HP08NC

For more information on any Conair Products call: 1-855-926-6247  
Visit our Web site: [www.conaircanada.ca](http://www.conaircanada.ca)  
Consumer Call Centre E-mail: [Consumer\\_Canada@Conair.com](mailto:Consumer_Canada@Conair.com)  
©2018 Conair Consumer Products ULC  
100 Conair Parkway, Woodbridge, ON L4H 0L2

IB-15905

18CN033692

BODY **BENEFITS**

by **CONAIR**



## HAVE A QUESTION?

*Please do not return this product to the retailer.*

**CALL US FIRST!**

Our customers service and product experts are ready to answer **ALL** your questions.

Please call our **TOLL-FREE** customer service number at :

**1-855-926-6247** Or visit us online at [www.conaircanada.ca](http://www.conaircanada.ca)

# Massaging Heating Pad

## Model HP08NC Instruction Booklet

For your safety and continued enjoyment of this products, always read the instruction book carefully before using.